# APPLICATION INSTRUCTIONS FOR “BIOREE” BIOLOGICAL SAMPLE STAINING KIT IN PREPARATION FOR SCANNING ELECTRON MICROSCOPY EXAMINATION

## 1. PURPOSE

1.1. “BioREE” reagent kit is designed for the preparation of biological samples for scanning electron microscopy examination.

## 2. KIT SPECIFICATIONS AND FUNCTIONING PRINCIPLE

2.1. Kit contents:

①	“1 – liquid for initial rinsing” (rinsing solution, labelled red) – 1 vial, 5 mL
②	“2 – contrasting by $\text{NdCl}_3+\text{NaCl}$ ” (first contrast agent solution, labelled yellow) – 1 vial, 2 mL
③	“3 – liquid for final rinsing” (final rinsing solution, labelled green) – 1 vial, 5 mL

2.2. “BioREE” kit is designed for simultaneous preparation of two samples.  
NOTE: if necessary, the kit may be split into two independent parts.

2.3. Kit functioning principle:

The initial rinsing removes the components of growth mediums and the liquids of tissue ground substances sorbed on the surface of the sample. The following soaking in the rare earth element solution leads to elective accumulation of the REE on cell membranes. The rinsing removes surpluses of the contrasting agent.

## 3. SAFETY PRECAUTIONS

3.1. All the components of the kit are non-toxic in used concentrations.

## 4. EQUIPMENT AND MATERIALS NECESSARY FOR USING THE KIT:

- receptacles for sample steeping in solutions;
- forceps for carrying the samples;
- a dropper capable of picking up to 2 mL of liquid volume for carrying the solutions to the staining container;

LLC “GLAUCON”, INN: 7722708109, OGRN: 1107746061114  
Correspondence address:  
Office 534, 11A Rossolimo str., Moscow, Russia, 119021

According to the patent application  
RU 2672356 C1.

[www.bioree.ru](http://www.bioree.ru)

- shaking machine capable of succussion with 3–4 mm amplitude at 4–8 Hz frequency at +37°C temperature (preferably);
- scanning electron microscope capable of detecting back-scattered electrons (BSE) is required for registering the results of staining.

## 5. PREPARING REAGENTS FOR THE ANALYSIS

5.1. Reagents must be soaked for 20 minutes at room temperature and then thoroughly mixed upturning every vial.

## 6. ANALYSIS PROCEDURE

6.1. Prepare appropriate amount of reagents depending on the number of samples being examined.

6.2. Put “1 – liquid for initial rinsing” into the first plate well. If the sample is adhered inside a container, it is possible to perform sample preparation by sequentially changing solutions in the container.

6.3. Carefully rinse the sample for 1–2 minutes.

6.4. Move the sample into the second plate well (container) and fill it with “2 – contrasting by  $\text{NdCl}_3 + \text{NaCl}$ ” agent.

6.5. Incubate the sample for 30 minutes at 37°C, preferably with succussion.

6.6. Move the sample into the third plate well (container) and fill it with “3 – liquid for final rinsing” solution.

6.7. Carefully rinse the sample.

6.8. Remove the surpluses of moisture.

6.9. Place the sample on the microscope stage. If examination cannot be carried out immediately, transportation and short-term storage of dried at RT samples are acceptable.

## 7. KIT STORAGE AND APPLICATION CONDITIONS

7.1. “BioREE-B” kit must be stored at +2° to +8°C temperature. Storage period is 2 months. Transportation and short-term storage at RT are acceptable.

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